

ASSEMBLY INSTRUCTION OA0-70A GCX Attachment for Monitor **FUKUDA DENSHI CO.,LTD.** 3-39-4, Hongo, Bunkyo-ku, Tokyo, Japan Tel:+81-3-5684-1455 Fax:+81-3-3814-1222



Copyright © 2014 by Fukuda Denshi Co., Ltd.
No part of this document may be copied or transmitted in any form without the prior written permission of Fukuda Denshi Co., Ltd.

- Thank you for purchasing this product.
- Before using this product, read through this assembly instruction and also refer to the DS-8200 operation manual.
- After reading, keep this assembly instruction near the unit for future reference.

CE This device bears the CE label in accordance with the provisions of Medical Device Directive 93/42/EEC.

EC REP Fukuda Denshi UK
Unit 7, Genesis Business Park, Albert Drive, Woking, Surrey GU21 5RW, United Kingdom

This product is used for mounting the Base Unit (BS-8210) of the DS-8200 monitor or the Countertop for DS-8200 (OA0-68A) on the GCX arm. Follow this assembly instruction, and use this product correctly.

Safety Precautions

Make sure to follow the precautions indicated below, as these are important messages related to safety. The warning displayed and symbols used in this instruction are shown below.

CAUTION	Failure to follow this message may cause injury or failure of the product.
NOTE	A note is not related to product safety, but provides information about the correct use and operating procedures to prevent incorrect operation and malfunction of the equipment.

NOTE


For quality improvement, specifications are subject to change without prior notice.

Contents

Refer to the following list and verify that all items are included.

OA0-70A x1	Screw M4X10 x4

Example



CAUTION

- This product is used when mounting the Base Unit (BS-8210) of the DS-8200 monitor or the Countertop for DS-8200 (OA0-68A) on the GCX arm. Do not use this product for any other purpose.
- Check that all the screws are tightened correctly after assembly and at periodic inspections.

Assembly Procedure

This procedure explains how to attach OA0-70A to BS-8210 or OA0-68A. The process requires the following tool.


- Phillips screwdriver #2

NOTE


- Attach this product to BS-8210 or OA0-68A on a level surface following the procedure below.
- A Phillips screwdriver is required to complete the installation.
- Make sure that the screws are securely tightened.

< Attached to BS-8210 >

1. Fix the OA0-70A on the bottom of the Base Unit (BS-8210) using the screws (4 locations).



2. Slide the OA0-70A on the top panel of the GCX arm and then secure the OA0-70A with the locking screws. (encircled in the following figure) ^{*1)}




^{*1)} For details, refer to the corresponding GCX operation manual.


< Attached to OA0-68A >

1. Fix the OA0-70A on the bottom of the Countertop for DS-8200 (OA0-68A) using the screws (4 locations).

* For details on how to attach the Base Unit (BS-8210) to the Countertop for DS-8200 (OA0-68A), refer to the OA0-68A Assembly Instruction.



2. Slide the OA0-70A on the top panel of the GCX arm and then secure the OA0-70A with the locking screws. (encircled in the following figure) ^{*1)}



^{*1)} For details, refer to the corresponding GCX operation manual.

Maintenance / Inspection

To prevent the equipment from falling off, verify that the screws are securely tightened once a year.